OT Glass II LC

Light-Cure Glass Ionomer Filling Cement Capsule For Professional Dental Use Only

Product Description

OT Glass II LC is a resin-modified glass ionomer restorative cement supplied in pre-dosed capsules for consistent performance and easy application. It combines the fluoride release and chemical adhesion of conventional glass ionomers with the strength, polishability, and esthetics of resin technology. The dual-cure system (light and chemical set) ensures complete polymerization even in deeper areas. Ideal for restorations, liners, and bases, OT Glass II LC provides durable sealing, long-term fluoride protection, and excellent marginal integrity in all clinical conditions.

Composition

Fluoro-aluminosilicate glass powder
Polyacrylic acid modified with HEMA
Bis-GMA and UDMA resins
Photoinitiators (camphorquinone, EDMAB)
Tartaric acid
Water
Pigments and stabilizers

Indications for Use

Class I, III, and V restorations.
Cervical lesions and root caries.
Small Class II restorations and pediatric restorations.
Core build-ups and sandwich technique under composites.
Base or liner under direct restorative materials.

Contraindications

Known hypersensitivity to methacrylates or glass ionomer components.
Pulp exposure or direct contact with pulp tissue.
Use in areas where complete moisture control is not possible.

Warnings

Avoid contact with eyes, skin, and soft tissues. If contact occurs, rinse immediately with water.
Uncured material may cause irritation or allergic reaction.
Use rubber dam isolation where possible.
For dental professional use only.

Precautions

Do not expose capsules to strong light before use.
Use the capsule immediately after activation and mixing.
Ensure cavity surfaces are clean and lightly moist — not desiccated.
Do not use eugenol-containing materials; they may interfere with polymerization.
Wear gloves and avoid prolonged contact with uncured material.

Properties (Typical)

Setting mechanism: Dual-cure (light + chemical)
Working time: 2 min at 23 °C
Setting time: 3–4 min after placement
Compressive strength: > 120 MPa
Flexural strength: > 40 MPa
Radiopacity: ≥ 3 mm Al
Fluoride release: Sustained and long-term
Film thickness: ≤ 25 µm
Shade: Universal tooth color

Directions for Use

Shake capsule to loosen powder.
Activate capsule and press plunger fully to break the membrane.
Insert capsule into the mixer and mix for 8 seconds (4000–5000 rpm).

Place capsule in the applicator gun and dispense directly into the cavity.
Sculpt and contour as required.
Light-cure for 20 seconds from each side using a curing light (≥1000 mW/cm²).
Finish after complete setting.

Storage and Handling

Store at 10–25 °C in a dry, cool place.
Do not freeze or expose to direct sunlight.
Use immediately after activation.
Shelf life: 24 months from manufacturing date (see package).

Adverse Reactions

Rare allergic responses to methacrylate or glass ionomer components may occur.
Discontinue use and seek medical advice if irritation, redness, or swelling appears.

Packaging

- Fifty (50) × pre-dosed capsules
- Instruction leaflet (this IFU)

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