OTFR Flow Composite EN

Flowable Fiber-Reinforced Composite For Professional Dental Use Only

Product Description

OTFR Flow Composite is an advanced **fiber-reinforced flowable restorative composite** that combines the adaptability of a flowable resin with the strength of randomly dispersed short glass fibers. The integrated fibers act as crack stoppers and stress absorbers, significantly enhancing toughness, fracture resistance, and longevity compared to conventional flowable composites.

With its optimized viscosity, OTFR Flow Composite flows easily into cavity walls and irregularities while forming a reinforced substructure under conventional composites. Designed for stress-bearing restorations, liners, and minimally invasive procedures, it ensures long-term reinforcement and protection in both posterior and cervical restorations.

Composition

Bis-GMA (Bisphenol A-glycidyl methacrylate)
UDMA (Urethane dimethacrylate)
TEGDMA (Triethylene glycol dimethacrylate)
Short glass fibers (E-glass, silanated) ≈ 25 wt%
Barium–alumino–borosilicate glass filler ≈ 50 wt%
Silica nano-filler ≈ 5 wt%
Photoinitiators (Camphorquinone, EDMAB)
Pigments and stabilizers

Indications for Use

Base or liner under direct composite restorations.
Reinforcement layer in stress-bearing restorations (Class I–II).
Minimally invasive posterior and cervical restorations.
Core build-ups and repair of fractured or defective restorations.
Substructure layer beneath universal or nanohybrid composites.

Contraindications

Known allergy or hypersensitivity to methacrylate resins or glass fibers.
Patients with a known sensitivity to (meth)acrylates.

Warnings

Avoid contact with eyes, skin, or mucosa. In case of contact, rinse immediately with water.
Uncured resin may cause irritation or sensitization.
Use rubber dam isolation when possible.
For dental professional use only.

Precautions

Ensure the cavity is clean, disinfected, and free from debris.
Do not use eugenol-containing materials prior to application.
Apply in increments ≤ 2 mm to ensure full polymerization.
Use a curing light with intensity ≥ 1000 mW/cm².
Do not polish the fiber-reinforced layer directly; overlay

with a conventional composite for optimal surface finish and esthetics.

Replace syringe cap immediately after use.

Properties (Typical)

Filler + fiber content: ≈ 82 wt% (≈ 61 vol%)
Working time: Unlimited (until light exposure)
Curing time: 20 seconds (2 mm layer, ≥ 1000 mW/cm²)
Flexural strength: > 160 MPa
Fracture toughness: > 2.5 MPa·m^{1/2}
Compressive strength: > 350 MPa
Elastic modulus: 10–12 GPa
Shrinkage: < 2 %
Radiopacity: ≥ 3 mm Al
Shade: Universal tooth color (translucent)
Viscosity: Thixotropic — flows under pressure, stable after placement

Directions for Use

- Tooth Preparation** — Prepare the cavity using standard restorative techniques. Clean and dry gently; avoid overdrying dentin.
- Bonding** — Apply a compatible universal or self-etch adhesive following the manufacturer's instructions.
- Application** — Dispense OTFR Flow Composite directly into the cavity using a dispensing tip. Spread evenly using a probe or brush to ensure adaptation to all walls.
- Light Curing** — Light-cure each increment (≤ 2 mm) for 20 seconds with an LED curing light (≥ 1000 mW/cm²).
- Capping Layer** — For optimal surface finish and polishability, overlay the reinforced layer with a conventional nanohybrid or universal composite.
- Finishing** — Finish and polish the final restoration using standard composite polishing systems.

Storage and Handling

Store at 10–25 °C in a cool, dry place.
Avoid exposure to direct sunlight or heat.
Do not freeze.
Keep syringe tightly capped when not in use.
Shelf life: 24 months from manufacturing date (see package).

Adverse Reactions

Rare allergic responses to methacrylate monomers or glass fiber components may occur.
Discontinue use and seek medical advice if redness, swelling, or irritation develops.

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